



Real-Time Excellence

see **MORE**
see **BETTER**
see **NOW**



Imagine a 3D Movement Analysis system that:

- Doesn't suffer from occlusions and never forgets which marker is which
- Is so precise, you don't need to filter data or show 'average' results; just show what you actually measured!
- Can be set up anywhere, not confined to a fixed laboratory
- Can be broken down and rebuilt in a different place in a new configuration – in ten minutes!
- Requires no calibration
- Shows your results immediately

...you're thinking about Codamotion



Codamotion CX1 System

It's easier

- No calibration
- No filtering
- No waiting for results

It's smaller

- Ultra wide viewing angle brings sensor unit closer to the measurement space (no wasted space)
- Utilises the latest miniaturised battery technology for markers

It's faster

- No post processing – calculated data analysis in seconds
- Customised marker systems for rapid, precise placement

It's complete

- Available as a complete kit including software, cabling, tripod and protective case

